


Cyanide A110 (w/In-line Dist. and UV Digest)


Example 1 - Range: 0.0005 - 0.2 mg/L



Cyanide A110-A01
Range: 0.0005 - 0.2 mg/L
MDL=0.0002 mg/L
Analysis rate: 30 samples/hour

NOTE: 0.0005 - 0.002 mg/L peaks are shown

Example 2 - Range: 0.005 - 1 mg/L



Cyanide A110-A00
Range: 0.005 - 1 mg/L
MDL=0.001 mg/L
Analysis rate: 40 samples/hour