


Cyanide A111 (Post Distillation)

Example 1 - Range: 0.0005 - 0.50 mg/L



Cyanide A111-A01


Range: 0.0005 - 0.50 mg/L

MDL=0.0001 mg/L

Analysis rate: 30 samples/hour

NOTE: 0.0005 - 0.002 mg/L peaks are shown

Example 2 - Range: 0.001 - 1.0 mg/L



Cyanide A111-A00

Range: 0.001 - 1.0 mg/L

MDL=0.0004 mg/L

Analysis rate: 55 samples/hour